


Competence CP3202S



Key features. The optimum for all weighing tasks..

- Excellent, easy-to-read backlit display with especially large digits (16.5 mm)
- Large weighing pan
- ISO/GLP-compliant logging of calibration, adjustment and weighing data using an optional data printer or a PC
- Four digital filter levels allow the balance to be optimally adapted to the particular ambient conditions at the place of installation.
- Very easy to operate; all keys feature positive click action; 2 large tare keys; clear and logically arranged keypad
- RS-232 serial interface port
- Integrated mechanical components that effectively protect the weigh cell – from overloading or impact against the sides of the weighing pan.
- Monolithic weigh cell technology – rugged, precise, reliable – year after year
- Rugged housing construction and use of housing materials that are highly resistant to chemicals – as a result, the balance can stand up to tough daily use
- Built-in application programs: counting, weighing in %, net-total formulation; animal weighing; mass unit conversion by toggling between 2 units



Data printer



SartoConnect

Readability	0.01 g
Weighing capacity	3,200 g
Pan size	190x204 mm
Response time (avg.)	1.5 s
Housing (WxDxH)	213x342x90 mm

Key accessories	Order no.
Data printer with date, time and statistics functions	YDP03-OCE
Dust cover	6960CP03
Data transfer software SartoConnect with RS-232C standard cable (1.5 m);	YSC01L
Calibration weight (1x2,000 g; F1)	YCW6238-00