

## ME5 and SE2 Ultra-Micro- and Microbalances



### Premium Balances Designed to Meet the Highest Requirements

- Outstanding weighing technology featuring exceptionally fast response times
- Round, motorized glass draft shield
- Draft shield operation with learning function activated by user-friendly draft shield keys
- Ideal for use as measuring and test equipment in quality management systems
- SQmin function for displaying the minimum sample quantity according to the United States Pharmacopeia
- New S.U.R.E. function for dynamic display of the uncertainty of measurement according to the DKD Calibration Certificate
- Fully automatic calibration|adjustment function, isoCAL
- Backlit graphical display with plain-text prompts and info on all settings
- Highly resistant membrane keypad overlay with "bubble-shaped" keys
- Complete package of applications is standard
- RS-232C port plus additional Sartorius printer port

Models	ME5	SE2
Readability	1 µg	0.1 µg
Weighing capacity	5.1 g	2.1 g
Repeatability	± ≤ 1 µg	± ≤ 0.25 µg
Max. linearity	± ≤ 4 µg	± ≤ 0.9 µg
Avg. response time	10 s	10 s
Pan size	Ø 30 mm	Ø 20 mm



Data printer



SartoConnect

Key accessories	Order No.
Data printer with date, time and statistics functions	YDP03-OCE
Foot switch for draft shield, tare and print functions	YPE01RC
Data transfer software SartoConnect with RS-232C standard cable (1.5 m);	YSC01L
Weight set for air density determination	YSS35-00
Ionizing blower	YIB01-ODR
Balance table	YWT01