



# QUALITY CONTROL SERVICES

LABORATORY EQUIPMENT • SALES • SERVICE • CALIBRATION • REPAIRS  
2340 SE 11<sup>TH</sup> Ave. Portland, Oregon 97214 • Box 14831 Portland, Oregon 97293  
(503) 236-2712 • FAX (503) 235-2535 • www.qc-services.com



Account Name: \_\_\_\_\_

Purchase Order or CC: \_\_\_\_\_

Contact Name: \_\_\_\_\_

Ship to Address: \_\_\_\_\_

Phone: \_\_\_\_\_

Dept/Bldg/Rm: \_\_\_\_\_

E-mail: \_\_\_\_\_

City, State, Zip: \_\_\_\_\_

PDF Calibration Certificates:  Yes (Certificates will be e-mailed to the address above. Hard copies will not be included with pipettes.)

Calibration Interval:  3 Months  6 Months  9 Months  12 Months  Other (specify): \_\_\_\_\_

Pipettes Externally Cleaned:  DNA Away (Add \$15)  RNase Away (Add \$15)  DNA & RNase Away (Add \$25)  
(Price is per pipette)

Expedited Turn-Around-Time:  24 Hrs. (Add \$15 per pipette - Contact QCS for advanced approval prior to shipping)  
(Not to exceed 15 pipettes)

Pipettes Decontaminated:  Yes By: Signature First/Last Name: \_\_\_\_\_

2022 Pipette Calibration Service Offerings						
Quantity	Pipette	Plan Type	As Found	Inspection & PM*	As Left	List Price*
	Single Ch.	Certificate of Calibration (COC)	2 x 3 No Data	Yes	2 x 5 No Data	\$50.00
	8-Ch Multi	Certificate of Calibration (COC)	2 x 3 No Data	Yes	2 x 5 No Data	\$85.00
	12-Ch Multi	Certificate of Calibration (COC)	2 x 3 No Data	Yes	2 x 5 No Data	\$105.00
	Single Ch.	Certificate w/ Data (CWD)	2 X 3	Yes	2 x 5	\$60.00
	8-Ch Multi	Certificate w/ Data (CWD)	2 X 3	Yes	2 x 5	\$95.00
	12-Ch Multi	Certificate w/ Data (CWD)	2 X 3	Yes	2 x 5	\$115.00
	Single Ch.	Certificate w/ Data & Unc. (A2LA)	3 X 3	Yes	3 x 10	\$70.00
	8-Ch Multi	Certificate w/ Data & Unc. (A2LA)	3 X 3	Yes	3 x 10	\$105.00
	12-Ch Multi	Certificate w/ Data & Unc. (A2LA)	3 X 3	Yes	3 x 10	\$125.00

\* Inspection & preventive maintenance (PM) include the inspection, functional tests, liquid leak check, cleaning, replacement, and lubrication of the liquid end components as needed with manufacturer parts and lubricant(s). Repairs and replacement of some parts may be at an additional charge with customer's approval. Quantity discounts may apply for 11+ pipettes sent in for calibration. Contact QCS for more details.

Customer Comments: \_\_\_\_\_